

Pharmacy Refrigerators HYC-68



Audible and Visual Alarms

Audible buzzer and visible flashing light



Five Alarm Functions

High/low temperature, power failure, door ajar, sensor error, low battery



Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

HYC-68

Technical Data

Cabinet Type	Upright, solid door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	35

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	84	68
Electrical Current(A)	0.65	1.1

Dimensions

Capacity(L/Cu.Ft)	68/2.4
Net/Gross Weight(approx)	38/40(kg) 83.8/88.2(lbs)
Interior Dimensions(W*D*H)	415*385*505(mm) 16.3*15.2*19.9(in)
Exterior Dimensions(W*D*H)	495*580*660(mm) 19.5*22.8*26.0(in)
Packing Dimensions(W*D*H)	540*640*710(mm) 21.3*25.2*28.0(in)
Container Load(20'/40'/40'H)	108/228/228



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

Accessories

Foot	Y
Shelves /Basket	2/1

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-68A



**Audible and
Visual Alarms**
Audible buzzer and visible
flashing light




**Five Alarm
Functions**
High/low temperature,
power failure, door ajar,
sensor error, low battery



**Air Cooling,
Auto Defrost**
Better cabinet temperature
uniformity, auto defrost

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

HYC-68A

Technical Data

Cabinet Type	Upright, glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	35

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	84	68
Electrical Current(A)	0.65	1.1

Dimensions

Capacity(L/Cu.Ft)	68/2.4
Net/Gross Weight(approx)	38/40(kg) 83.8/88.2(lbs)
Interior Dimensions(W*D*H)	415*385*505(mm) 16.3*15.2*19.9(in)
Exterior Dimensions(W*D*H)	495*580*660(mm) 19.5*22.8*26.0(in)
Packing Dimensions(W*D*H)	540*640*710(mm) 21.3*25.2*28.0(in)
Container Load(20'/40'/40'H)	108/228/228



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

Accessories

Foot	Y
Shelves /Basket	2/1

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1

Haier Road, Qingdao, 266101, P.R. China

Tel: +86-0532-88936011 / 5955

Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-118



**Audible and
Visual Alarms**
Audible buzzer and visible
flashing light



**Five Alarm
Functions**
High/low temperature,
power failure, door ajar,
sensor error, low battery



**Air Cooling,
Auto Defrost**
Better cabinet temperature
uniformity, auto defrost

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

HYC-118

Technical Data

Cabinet Type	Upright, solid door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	32

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	135	75
Electrical Current(A)	1.03	1.11

Dimensions

Capacity(L/Cu.Ft)	118/4.2
Net/Gross Weight(approx)	41/46(kg) 90.4/101.4(lbs)
Interior Dimensions(W*D*H)	515*415*630(mm) 20.3*16.3*24.8(in)
Exterior Dimensions(W*D*H)	597*635*835(mm) 23.5*25.0*32.9(in)
Packing Dimensions(W*D*H)	680*690*910(mm) 26.8*27.2*35.8(in)
Container Load(20'/40'/40'H)	48/102/102



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	3/1

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Haier Biomedical

Intelligent Protection of Life Science

Pharmacy Refrigerators HYC-118A



**Air Cooling,
Auto Defrost**
Better cabinet temperature
uniformity, auto defrost



**Five Alarm
Functions**
High/low temperature,
power failure, door ajar,
sensor error, low battery



**Audible and
Visual Alarms**
Audible buzzer and visible
flashing light

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

HYC-118A

Technical Data

Cabinet Type	Upright, glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	32

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	165	140
Electrical Current(A)	1.16	1.46

Dimensions

Capacity(L/Cu.Ft)	118/4.2
Net/Gross Weight(approx)	46/51(kg) 101.4/112.4(lbs)
Interior Dimensions(W*D*H)	515*415*630(mm) 20.3*16.3*24.8(in)
Exterior Dimensions(W*D*H)	597*635*835(mm) 23.5*25.0*32.9(in)
Packing Dimensions(W*D*H)	680*690*910(mm) 26.8*27.2*35.8(in)
Container Load(20'/40'/40'H)	48/102/102



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	3/1

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-260



**Air Cooling,
Auto Defrost**
Better cabinet temperature
uniformity, auto defrost



**Five Alarm
Functions**
High/low temperature,
power failure, door ajar,
sensor error, low battery



Adjustable Shelves
The shelves are adjustable to
any height to fulfill different
requirements

Ergonomic Design



High Accuracy
Temperature Control



Interior
Fluorescent
Light



Lockable
Design

Specifications

HYC-260

Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	55

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	208~230/60
Power(W)	300	300
Electrical Current(A)	2.1	2.1

Dimensions

Capacity(L/Cu.Ft)	260/9.2
Net/Gross Weight(approx)	88/100(kg) 194.0/220.5(lbs)
Interior Dimensions(W*D*H)	550*460*1065(mm) 21.7*18.1*41.9(in)
Exterior Dimensions(W*D*H)	620*655*1720(mm) 24.4*25.8*67.7(in)
Packing Dimensions(W*D*H)	695*720*1795(mm) 27.4*28.3*70.7(in)
Container Load(20'/40'/40'H)	24/48/48



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Shelves /Basket	4/-

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-290



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Air Cooling, Auto Defrost
Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves
The shelves are adjustable to any height to
fulfill different requirements



USB Interface
Data of temperature history can be downloaded
via the built-in USB port and optional printer design



Heated Double Pane Glass Door
Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



**Self Closing Door with 90° Stay
Open Feature**
Energy conserving, maintains temperature
uniformity to protect valuable samples

Specifications

HYC-290

Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	43

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/60	115/60
Power(W)	340	460
Electrical Current(A)	2.2	5.0

Dimensions

Capacity(L/Cu.Ft)	290/10.2
Net/Gross Weight(approx)	105/127(kg) 231.5/280.0(lbs)
Interior Dimensions(W*D*H)	530*555*1080(mm) 20.9*21.9*42.5(in)
Exterior Dimensions(W*D*H)	665*710*1665(mm) 26.2*28.0*65.6(in)
Packing Dimensions(W*D*H)	735*800*1800(mm) 28.9*31.5*70.9(in)
Container Load(20'/40'/40'H)	21/42/42



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	5/-
USB Interface	Y

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-360



Five Alarm Functions

High/low temperature,
power failure, door ajar,
sensor error, low battery



Air Cooling, Auto Defrost

Better cabinet temperature
uniformity, auto defrost



Adjustable Shelves

The shelves are adjustable to
any height to fulfill different
requirements

Ergonomic Design



High Accuracy
Temperature Control



Interior
Fluorescent
Light



Lockable
Design

Specifications

HYC-360

Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	55

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	208~230/60	115/60
Power(W)	315	315	340
Electrical Current(A)	2.2	2.2	4.5

Dimensions

Capacity(L/Cu.Ft)	360/12.7
Net/Gross Weight(approx)	99/114(kg) 218.3/251.3(lbs)
Interior Dimensions(W*D*H)	550*460*1340(mm) 21.7*18.1*52.8(in)
Exterior Dimensions(W*D*H)	620*655*1995(mm) 24.4*25.8*78.5(in)
Packing Dimensions(W*D*H)	695*720*2080(mm) 27.4*28.3*81.9(in)
Container Load(20'/40'/40'H)	24/48/48



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Shelves /Basket	5/-

Others

Certificate	CE	/
-------------	----	---

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-360



Five Alarm Functions

High/low temperature,
power failure, door ajar,
sensor error, low battery



Air Cooling, Auto Defrost

Better cabinet temperature
uniformity, auto defrost



Adjustable Shelves

The shelves are adjustable to
any height to fulfill different
requirements

Ergonomic Design



High Accuracy
Temperature Control



Interior
Fluorescent
Light



Lockable
Design

Specifications

HYC-360

Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	55

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	208~230/60	115/60
Power(W)	315	315	340
Electrical Current(A)	2.2	2.2	4.5

Dimensions

Capacity(L/Cu.Ft)	360/12.7
Net/Gross Weight(approx)	99/114(kg) 218.3/251.3(lbs)
Interior Dimensions(W*D*H)	550*460*1340(mm) 21.7*18.1*52.8(in)
Exterior Dimensions(W*D*H)	620*655*1995(mm) 24.4*25.8*78.5(in)
Packing Dimensions(W*D*H)	695*720*2080(mm) 27.4*28.3*81.9(in)
Container Load(20'/40'/40'H)	24/48/48



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Shelves /Basket	5/-

Others

Certificate	CE	/
-------------	----	---

Haier Biomedical UK

Address: 2 Birch House, Brambleside, Bellbrook
Industrial Estate, Uckfield, East Sussex,
TN22 1QQ, United Kingdom
TEL: +44(0) 1825 744 696 (Head Office)
www.haierbiomedical.co.uk



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Scan 2D barcode with your mobile device to learn more information about Haier Biomedical

Pharmacy Refrigerators HYC-390



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Air Cooling, Auto Defrost
Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves
The shelves are adjustable to any height to
fulfill different requirements



USB Interface
Data of temperature history can be downloaded
via the built-in USB port and optional printer design



Heated Double Pane Glass Door
Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



**Self Closing Door with 90° Stay
Open Feature**
Energy conserving, maintains temperature
uniformity to protect valuable samples

Specifications

HYC-390

Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	43

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/60	115/60
Power(W)	380	500
Electrical Current(A)	2.4	5.5

Dimensions

Capacity(L/Cu.Ft)	390/13.8
Net/Gross Weight(approx)	116/139(kg) 255.7/306.4(lbs)
Interior Dimensions(W*D*H)	530*555*1380(mm) 20.9*21.9*54.3(in)
Exterior Dimensions(W*D*H)	665*710*1965(mm) 26.2*28.0*77.4(in)
Packing Dimensions(W*D*H)	735*800*2100(mm) 28.9*31.5*82.7(in)
Container Load(20'/40'/40'H)	21/42/42



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	7/-
USB Interface	Optional

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-390F



Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements



USB Interface

Data of any height to fulfill different requirements via the built-in USB port and optional printer design



Self Closing Door with 90° Stay Open Feature

Energy conserving, maintains temperature uniformity to protect valuable samples

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

HYC-390F

Technical Data

Cabinet Type	Upright, single solid door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	43

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/60
Power(W)	300
Electrical Current(A)	2.4

Dimensions

Capacity(L/Cu.Ft)	390/13.8
Net/Gross Weight(approx)	106/129(kg) 233.7/284.4(lbs)
Interior Dimensions(W*D*H)	530*555*1380(mm) 20.9*21.9*54.3(in)
Exterior Dimensions(W*D*H)	665*710*1965(mm) 26.2*28.0*77.4(in)
Packing Dimensions(W*D*H)	735*800*2100(mm) 28.9*31.5*82.7(in)
Container Load(20'/40'/40'H)	21/42/42



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	7/-
USB Interface	Y

Others

Certificate	CE,UL
-------------	-------

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-610



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Heated Double Pane Glass Door

Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



Self Closing Door with 90° Stay Open Feature

Energy conserving, maintains temperature
uniformity to protect valuable samples

Specifications

HYC-610

Technical Data

Cabinet Type	Upright, single door
Climate Class	ST
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	43

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220~240/50/60	115/60
Power(W)	550	560
Electrical Current(A)	3.5	7.5

Dimensions

Capacity(L/Cu.Ft)	610/21.6
Net/Gross Weight(approx)	204/227(kg) 449.7/500.4(lbs)
Interior Dimensions(W*D*H)	680*640*1400(mm) 26.8*25.2*55.1(in)
Exterior Dimensions(W*D*H)	780*840*1945(mm) 30.7*33.1*76.6(in)
Packing Dimensions(W*D*H)	865*945*2090(mm) 34.1*37.2*82.3(in)
Container Load(20'/40'/40'H)	12/26/26



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	6/-
USB Interface	Optional
Temp Recorder	Optional

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-940



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Heated Double Pane Glass Door

Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



Self Closing Door

Energy conserving, maintains temperature
uniformity to protect valuable samples

Specifications

HYC-940

Technical Data

Cabinet Type	Upright, double glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	45

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	208~230/60	115/60
Power(W)	850	770	850
Electrical Current(A)	4.5	4.5	11

Dimensions

Capacity(L/Cu.Ft)	890/31.4
Net/Gross Weight(approx)	227/279(kg) 500.4/615.1(lbs)
Interior Dimensions(W*D*H)	1030*590*1425(mm) 40.6*23.2*56.1(in)
Exterior Dimensions(W*D*H)	1130*755*1980(mm) 44.5*29.7*78.0(in)
Packing Dimensions(W*D*H)	1250*850*2130(mm) 49.2*33.5*83.9(in)
Container Load(20'/40'/40'H)	8/18/18



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	12/-
USB Interface	Optional

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-940F



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Self Closing Door

Energy conserving, maintains temperature
uniformity to protect valuable samples

Specifications

HYC-940F

Technical Data

Cabinet Type	Upright, double solid door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	45

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
Power(W)	600
Electrical Current(A)	3

Dimensions

Capacity(L/Cu.Ft)	890/31.4
Net/Gross Weight(approx)	207/259(kg) 456.4/571.0(lbs)
Interior Dimensions(W*D*H)	1030*590*1425(mm) 40.6*23.2*56.1(in)
Exterior Dimensions(W*D*H)	1130*755*1980(mm) 44.5*29.7*78.0(in)
Packing Dimensions(W*D*H)	1250*850*2130(mm) 49.2*33.5*83.9(in)
Container Load(20'/40'/40'H)	8/18/18



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	12/-
USB Interface	Optional

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Pharmacy Refrigerators HYC-1378



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Special Designed Slot for Sign Card

Adjustable shelves, the slot design be clear,
store items is easy to distinguish



Heated Double Pane Glass Door

Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



Self Closing Door

Energy conserving , maintains temperature
uniformity to protect valuable samples

Specifications

HYC-1378

Technical Data

Cabinet Type	Upright, double door		
Climate Class	N		
Cooling Type	Forced air cooling		
Defrost Mode	Auto		
Refrigerant	CFC-Free		
Noise((dB(A))	45	45	48

Performance

Temp Range(°C)	2~8		
----------------	-----	--	--

Control

Controller	Microprocessor		
Display	LED		

Electrical Data

Power Supply(V/Hz)	220-240/50	208~230/60	115/60
Power(W)	800	800	950
Electrical Current(A)	4.2	4.2	11

Dimensions

Capacity(L/Cu.Ft)	1378/48.7		
Net/Gross Weight(approx)	310/360(kg) 683.4/793.7(lbs)		
Interior Dimensions(W*D*H)	1320*700*1500(mm) 52.0*27.6*59.1(in)		
Exterior Dimensions(W*D*H)	1440*925*1980(mm) 56.7*36.4*78.0(in)		
Packing Dimensions(W*D*H)	1570*1010*2120(mm) 61.8*39.8*83.5(in)		
Container Load(20'/40'/40'H)	6/12/12		



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	12/-
USB Interface	Y
Temp Recorder	Optional

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.

